

 Certified Refurbished

## The smart choice for reliable anesthesia

The Narkomed 2B is one of Dräger's longest running and most reliable anesthesia machines. It comes standard with O<sub>2</sub> and N<sub>2</sub>O and holds three vaporizers. The Narkomed 2B also includes F-Style ascending bellows with oxygen, volume and pressure monitoring.

### A Quick Look

- » 19.1 and 19.3 vaporizers available
- » ORM standard; ORMC optional
- » Ventilator options: Pneumatic or electronic, ascending or descending bellows
- » Oxygen monitor and CO<sub>2</sub> absorber standard
- » AVE Electronic ventilator
- » Absorber and scavenger
- » Oxygen analyzer
- » Breathing pressure monitor
- » Respiratory volume
- » Centralized alarm and data
- » Oxygen ratio controller
- » Pressure limit controller



*Drawer pack configurations may vary from picture.*

### Technical Specifications

Maximum Height .....	68" (173 cm)
Maximum Width .....	40" (102 cm)
Maximum Depth .....	25" (63.5 cm)
Table-top Height .....	34" (86 cm)
Table-top Area .....	271 in. sq. (1,749 cm <sup>2</sup> )
Shelf Area .....	200 in. sq. (1,290 cm <sup>2</sup> )
Drawer Dimensions	
- Large .....	14"x21"x7" (36 x 52 x 18 cm)
- Small .....	14"x17"x3" (36 x 43 x 8 cm)
Nominal Weight .....	325 lb. (147 kg)