


Holds three vaporizers for multiple agents

- » Oxygen, volume and pressure monitoring on a removable flat plate display
- » O₂ and N₂O standard; air optional
- » ORC (Hypoxic Guard)
- » AVE ventilator with ascending bellows
- » Compatible with 19.1 or 19.3 vaporizers (optional)
- » E-Cylinder and central gas connections
- » CO₂ absorber standard



Specifications

- » Weight: Approximately 400 lbs.
- » Maximum dimensions: 40" x 68" x 25" (l x h x w)
- » Top shelf: 32" x 15" (l x w)
- » Table top: 271 square inches



Drawer pack configurations may vary from picture