

 **DRE Milano E20**  
Power Procedure Chair



Engineering with **Elegance**

# Power Procedure Chair

Naturally adapts to almost every physician, patient, procedure and environment

The DRE Milano E20 is an excellent alternative to traditional procedure chairs that offers innovative features and modern styling at a competitive price. The Milano E20 is comfortable for both the operator and patient: it is wider than traditional procedure chairs and offers easy access to people with mobility difficulties, thanks to the chair's adjustable height. Plus, its physician-configured adjustments provide excellent visibility within the operative field.

***Designed for comfort, the Milano E20 can be used for lengthy office-based procedures and is perfect for Plastic Surgery, Dermatology, Oral & Maxillofacial, ENT / Otolaryngology, Ophthalmology, OB/GYN and Diagnostic / Ultrasound Imaging.***



Multiple electric movements and adjustable height offer easy access for patients



## Standard Accessories



Handset Remote Control

Foot Pedal Control



Moveable Arm Rests



Paper Roll Holder




Seat Rails



## Optional Accessories

Instrument Tray



Castors



Articulated Knee Crutch





## Specifications:

- Maximum Capacity Load: 440 lbs. (200 kg)
- Angle of inclination: Trendelenburg: 30°, Back rest: 80°, Leg rest: 90°, Rotation (both sides): 90°
- Height: Minimum - 20" (52 cm), Maximum - 32" (82 cm)
- Dimensions: 67" x 30" x 34" (170 x 74 x 86 cm)
- Mattress dimensions: Length - 66" (168 cm), with extended leg rest - 73" (185 cm); Width - 23" (59 cm)
- Couch Weight: 168 lbs. (75 kg)
- Voltage: 100-110 - 230V
- Latex-free, non-toxic and non-carcinogenic materials
- Fire-resistant upholstery: UNI 9175 class 1M

