

DRE Milano T50 Power Procedure Table



Engineering with **Elegance**

Power Procedure Table

Naturally adapts to almost every physician, patient, procedure and environment

The DRE Milano T50 is an electric multi-function table with four actuators capable of adjusting height, backrest, legrest and trendelenburg with high precision and ensuring the patient an easy access in and out of the chair. It also features an automatic reset switch that returns the table to its starting position.

Designed for comfort, the Milano T50 can be used for lengthy office-based procedures and is perfect for Plastic Surgery, Dermatology, Oral & Maxillofacial, ENT / Otolaryngology, Ophthalmology, OB/GYN and Diagnostic / Ultrasound Imaging.





Optional Accessories



Articulating Arm Support

Narrow Headrest



Articulating Knee Crutch



Instrument Tray



Medical Armrest



I.V. Pole



Curtain Holder



Casters



Specifications:

- Maximum Capacity Load: 440 lbs. (200 kg)
- Height: Minimum - 24" (62 cm), Maximum - 36" (92 cm)
- Angle of inclination: Trendelenburg: 25°, Back rest: 75°, Leg rest: 65°
- Latex-free, non-toxic and non-carcinogenic materials
- Dimensions: 74" x 24" x 34" (190 x 62 x 86 cm)
- Weight: 196 lbs. (89 kg)
- Voltage: 110 and 220 V cc
- Warranty: Table – 2 year warranty; Mattress – 5 year warranty

