

## Comfortable, Functional and Economical

This versatile table is a perfect option for family practice, doctor's offices, clinics or wherever a cost effective treatment or exam table is needed. The DRE Patient Exam Table maximizes convenient features including built-in stirrups, leg rest, and step stool that can be easily pulled out or tucked away.

### Features:

- Durable vinyl upholstery over two-inch firm padding.
- Full 27-inch width for patient comfort.
- Air-spring adjustable backrest.
- Coated, two-position adjustable stirrups.
- Padded, pull-out leg rest on steel track
- Easy care all laminate base.
- Chip resistant edges.
- One large side storage compartment with shelf.
- Storage compartment is accessible from both sides.
- Two Euro style drawers with all steel slides.
- 400 lbs. load capacity under normal use.
- Plywood and laminate panels are California CARB I compliant.
- Top sections are easily replaceable in the field.
- Adjustable feet for leveling.
- Optional — convenient, pull-out step stool constructed with welded steel.
- Optional — pull-out basin and paper cutter.

### Dimensions:

- Exam Table: 58 in (72 in extended) L x 28.5 in W x 31 in H
- Step Stool (Optional): 16.875 in L x 19.5 in W x 7.75 in H



88-01-2016-09