

Reliable, Economical, Effective

The DRE SLS 6000 Surgery Light halogen light has been proven effective through years of use in operating rooms and features adjustable illumination, a high color temperature, and multiple mounting options to make it ideal for surgical procedures.

Features:

- OSRAM Halogen Bulbs
- Polished panel
- Removable handle
- Automatic switch within 0.2 seconds

Specifications:

- Illumination: 40,000-80,000 Lux (adjustable)
- Diameter of light spot: 120-300mm
- Depth of illumination: 1200mm
- Color temperature: 4500K±500K
- Increase of temperature: ≤10°C
- Power: 110W
- Lamp holder motion radius: 1800mm
- Multi-mirror: 2600 pieces
- CRI: ≥95%



Control Panel



Available as portable unit, single or dual- head ceiling mount, or single-head wall mount