

## Energy saving and eco-friendly, with an ultra-long service life

The DRE SLS 9000 Portable LED Light is a cost efficient portable light boasting a rated LED life of up to 30,000 to 50,000 hours.

The CAD designed multi-LED-lighting-beam produces lighting depth more than 55in., with up to 90,000 lux of central illumination. New alloy construction and a simple suspension system is easy to control and to be positioned accurately. This easily adjustable light system is suitable for pendant, wall-mounted and floor-standing lamps, meeting requirements of multiple operating rooms.



### Features:

- OSRAM led bulbs
- Adjustable focusing
- Diameter (in): 20
- Illumination (Lx):=90,000
- Color Temp (K): 3500-5000
- Color Rendering Index (%): =93
- Diameter of light spot (in): 5-14
- LED power (W): 55
- Illumination Depth (in): 28-59
- LED service life (H): 50,000



Easy Control Panel