

DRE Torino

Orthopedic Surgery Table

Equipment for the way *you* operate

The **DRE Torino** Orthopedic Surgery Table is the perfect choice to perform traction procedures. The Tractor is to be used in conjunction with the Torino EXL Surgical Table and work seamlessly in the Operating Room.



Technical Specifications

- » Height \leq 800mm; 32 in.
- » Lifting range \geq 350mm; 14 in.
- » Traction range \geq 190mm; 7 in.
- » Traction angle, rotation(up & down) 120°
- » Horizontal rotation 360°
- » Horizontal rotation of body frame \geq 180°
- » Extendable adjustment (back & forth) \geq 400mm; 16 in.
- » Lifting adjustment range of side frame \geq 200mm; 8 in.

