

## The **automatic** DRE Om2100 is perfect for cataract surgeries and post-op exams

- Ideal for precision procedures
- Used by ophthalmologists and surgeons in demanding environments
- Compact and easy-to-position
- Features motorized zoom and focus
- Comes standard with a new light source, fiber-optic cable, power supply and portable four-wheel base



### ► A Quick Look

- Motorized focus and zoom including foot pedal
- 65 mm interchangeable objective lens
- New fiber-optic illumination, including fiber-optic cable
- Zeiss-compatible upgradable components
- Inclined binoculars
- Mark II light source
- Mark II light tray
- Observation kit
- Beam splitter

Includes a **foot pedal** for hands-free focus and zoom