

Get flexibility and durability with DRE's manual ENT microscope

Featuring flexible, manual magnification, the DRE Em1000 is perfect for office settings. Its durable, easy-to-position stand lets you raise and lower the microscope without changing the horizontal position of the head, maintaining the focus of the microscope.



Choose between two types of bases

European-style tri-base (left) • Portable four-wheel base (right)



► A Quick Look

- Three-step manual magnification changer
- 48 mm interchangeable objective lens
- New fiber-optic illumination, including fiber-optic cable
- Zeiss-compatible upgradable components
- New 150-watt halogen light source