

DRE's MRI-compatible anesthesia machine is designed to meet all your anesthesia requirements

- » Three gas machine (O₂ • N₂O • Air).
- » A smart alternative to refurbished machines.
- » Advanced easy-to-use features for maximum patient safety.
- » Holds two vaporizers.
- » Large work area with a writing table and footrest.



A Quick Look

- » MRI-compatible.
- » Ventilator option (not included): Nuffield MRI-compatible anesthesia ventilator.
- » Selectatec®-compatible vaporizer mounting system.
- » Holds up to two vaporizers.
- » Specifically designed for low flow.
- » Left or right-handed layout.
- » Contemporary ergonomic design.
- » Advanced easy-to-use features for maximum patient safety.
- » Mechanical anti-hypoxic device and air/N₂O interlock.
- » Only two preventive maintenances required each year.

Physical Dimensions

- » Size (h x w x d): 54" (137 cm) x 26" (66 cm) x 26" (66 cm).
- » Weight (with all options): 200 lbs. (75 kg).

Power

- » Power Requirements: 110/220 VAC.
- » Auxiliary Electric Outlets: (4) 13 amp.

Display

- » Top Shelf: 24" (61 cm) x 16" (40 cm).
- » Work Surface: 24" (61 cm) x 14.5" (37 cm).
- » Drawers: 6" (15 cm) x 15" (38 cm) x 17" (43 cm), maximum of three.
- » Writing Tablet: A4 size.

Available with the Nuffield Anesthesia Ventilator

- » MRI-compatible, pneumatically driven time-cycled ventilator.
- » Pre-set volume and flow rate for adult or pediatric patients.
- » Powered by medical air or oxygen at 340 kPa - 410 kPa.