

Versatility for a more efficient Endoscopy Department

Designed for every stage of endoscopy — carries patients from pre-op to procedure to recovery.

Offers **C-arm compatibility** for imaging and the benefits of a **fully-featured surgical table**.

Improves safety for your staff — eliminates the need for patient transfers, reducing the risk of injury to your staff and the patient!

Ergonomic design — narrow stretcher width and high clearance base provide clinicians with comfortable access to patients.

Stretcher can be positioned as a chair for patient recovery.

Standard surgical side rail for attachment of accessories.


Easy to steer and control.

Full radiolucent top.



Lowers to 24" for universal patient access and **raises to 40"** for physician comfort during procedures.

Patients can remain on the Versailles from admission to discharge — even when imaging or surgery is required!



The singular endoscopy solution



Patient Transport
Moves patients through the endoscopy unit.



Chair Positioning
Provides patient comfort during recovery.



Adjustable Height
Low: 24"
High: 40"



Four Independent Safety Rails
You can keep the shoulder rail in place for patient support, and lower the hip rail for access to the patient.



Remote Control
Powered version only.



C-arm Access
Powered version: 100%
Manual version: 75%



Surgery Table Positioning Built-in



Removable Head Section



Push button rotating cot sides



Accessories available for alternative positioning



Removable ultra lightweight leg section



Mattress conforms to patient for increased comfort

Specifications	Powered (V100)	Manual (M75)	Additional Features
Height (min)	23.8" (605 mm)	22.6" (575 mm)	Lateral Tilt*, Head End Traverse*, Foot End Traverse*, K8 Pressure Care Mattress, Removable Ultra-Lightweight Leg Section, Removable and Rotating Safety Rails, Push Button Release/Self Locking Controls, 5th Wheel Steering Mechanism, Advanced Four-Caster Braking System, E Size Oxygen Cylinder Trough (Side Loading), Transfusion Pole with Mounting Position, and 150 mm Casters
Height (max)	39.5" (1,005 mm)	36.8" (935 mm)	
Length	82.7" (2,100 mm)	82.7" (2,100 mm)	
Width without Safety Rail	25.8" (655 mm)	25.8" (655 mm)	
Width with Safety Rail	30.3" (770 mm)	30.3" (770 mm)	
Mattress Length	80.3" (2,040 mm)	80.3" (2,040 mm)	
Mattress Width	23.6" (600 mm)	23.6" (600 mm)	
Mattress Depth	3" (75 mm)	3" (75 mm)	
Base Length	47" (1,190 mm)	41.5" (1,055 mm)	
Trendelenburg	20°	20°	
Reverse Trendelenburg	12°	N/A	
Lateral Tilt	+12°/-12°	N/A	
Backrest Articulation	0 - 65°	0 - 80°	
Head Section Articulation	+30°/-30°	+30°/-30°	
Head End Traverse	9" (230 mm)	N/A	
Foot End Traverse	9" (230 mm)	N/A	
C-arm Access	100%	75%	
Max. Patient Weight Limits	550 lbs (250 kg)	352 lbs (160 kg)	

Specifications subject to change without notice.

**Powered configuration only*