

Bright, clean, and even field of illumination



The **DRE Xavier C3 Portable LED** headlight provides a bright, clean and even field of illumination. With all of the features of the DRE Xavier C2, the C3 has the added intensity, longevity and brilliance that LED technology now provides.

The Xavier C3 Portable LED applies advanced High Beam TriLens technology for the result of unequaled clean, bright, adjustable and well-focused large spot of light.

The Xavier C3 Lithium Battery Pack with Intensity Control provides 7.5 hours of highest intensity light with up to 600 recharges.

You can now upgrade any of your older model DRE or BFW halogen headlights to LED with our conversion kit. Easily and affordably convert your light to LED with a new LED light harness, lithium ion battery pack and charger. As long as your headlight was manufactured after 1995 with the rectangular bulb housing, the conversion will work.



Features:

- 50,000 hour 3-watt LED
- 6,000° K spot (3,200° K LED optional)
- 4,500 foot candles at 10" working distance
- Variable intensity control
- Light spot size from .25" to 2" at a 10" working distance
- Battery Life of 7.5 to 100 hours, based on intensity
- Quick, 3-hour battery charge time
- Pinpoint spot adjustment
- Lightweight, comfortable headband (Options for firm nylon or fabric sport band)
- 100v-240v charger
- Autoclavable positioning bar provides sterile adjustability
- Pivot points allow for exceptional coaxial alignment



DRE Xavier C3 Portable LED

The full system shown here with LED headlight, battery pack and charger.