

Combination of classic aesthetics and modern technology

The GE Giraffe OmniBed combines two traditional infant care products in one. A full-featured incubator and a radiant warmer in one. Protects from heat loss while the bi-directional maintains thermal stability. Infant receives uniform heat regardless of position.

Features:

- Incubator with an Air Curtain to protect from heat loss
- Radiant warmer that provides consistent heat in any position
- A consistent and controlled thermal environment
- Improved access to patient
- A less stressful environment for the patient and family
- No need to transfer patient from incubator to infant warmer
- 360 degree rotating patient mattress
- Drop down removable doors and rotating bed provide easy access
- In-Bed scale



Certified Refurbished



Specifications:

- Dimensions
 - » Maximum height: Canopy down: 70 in (178 cm)
Canopy up: 94 in (239 cm)
 - » Minimum height: Canopy down: 58 in (147 cm)
Canopy up: 82 in (208 cm)
 - » Footprint: 45 in × 26 in (114 cm × 66 cm)
 - » Weight: 284 lbs. (129 kg)
 - » Floor to mattress height: 32 in to 44 in (81 cm to 111 cm)
 - » Mattress size: 26 in × 19 in (66 cm × 48 cm)
 - » Access door size: 9 in × 31 in (23 cm × 78 cm)
 - » Drawer size: 19 in × 20 in (48 cm × 50 cm)
 - » Drawer depth: 8 in (20 cm)
- Physical Characteristics
 - » Mattress tilt angle: 12°, continuously variable
 - » Tubing access ports: 8
 - » Servo humidity control range: 30-95%
 - » Relative humidity in 5% increments
 - » Alarm sound level: Adjustable audible levels
- In-Bed Scale Performance
 - » Accuracy: ±0.35 oz (10 g)
 - » Range: 0.66 lbs to 17.6 lbs (300 g to 8 kg)