

## Refurbished Hill-Rom GPS 880 Series stretchers provide smooth mobility and support weights up to 500 pounds.

Perfect for hospitals, surgery centers and teaching facilities!

### Specifications

- Overall length: 84"
- 8" casters
- Weight capacity: 500 lbs.
- Trendelenburg/Reverse Trendelenburg
- Dual-end brake and steer
- Sleep surface: 26" (w) x 78" (l)
- Height range: 23.5" to 28.5"
- Four IV rod locations
- One permanent IV rod
- Integrated storage bin
- Full perimeter bumper
- Stow-away side rails
- Radiolucent surface

### Options

(subject to availability)

- Oxygen bottle holder
- Pneumatic head section
- Pneumatic thigh gatch
- 5th wheel steering

### Refurbishing Procedures

- New paint and labels
- New rubber bumper
- New pedals
- Fully rebuilt hydraulic cylinder
- New gas dampeners
- Rebuilt casters with new components and wheels
- Full function tests

