

Hill-Rom Hybrid

For safety, quality and comfort.

Intelligent, well designed, clinically-effective and easy to use.

Our Process Includes

- Professionally refurbished medical beds to manufacturer's specifications
- JCAHO required documentation
- New paint & labels
- Fully function tested
- Basic bed function (Head, Knee, Hi-low)
- Orthopedic/Hard pan surface
- Reconditioned PC boards, 1 year warranty
- Fully rebuilt side rails w/new plastics and switches
- New head/foot boards, your choice of laminate colors, lifetime warranty
- New paint Light Neutral
- New motor couplings and capacitors
- Rebuilt casters with new components and wheels
- 1 year warranty (Domestic only)

Specifications

- 4-Motor bed
- Hi-Lo travel time: 21 seconds
- Sleep surface: 36"x 82"
- Overall length (w/out bumpers): 36" x 91"
- Bed frame clearance w/5" casters: 7.5"
- Low position: 17"
- High position: 30"
- Central brake and steer
- Automatic night light
- Instant CPR release

Options

- Radio/TV/ Lighting
- Bed exit
- Nurse call
- New PC boards, 5 year warranty
- Parallelogram (new style) side rails
- Zone Aire
- Dynamic Aire
- Mattress (built to OEM spec)
- Service Contracts for preventative maintenance
- 220-volt transformers

