

Hill-Rom VersaCare

For safety, quality and comfort.

Intelligent, well designed, clinically-effective and easy to use.

Our Process Includes

- Professionally refurbished medical beds to manufacturer's specifications
- JCAHO required documentation
- New paint & labels
- Fully function tested
- Reconditioned PC boards and side rail switches
- Reconditioned side rails
- Reconditioned head/foot boards
- New led plus trend indicat or
- Fresh hydraulic fluid
- Rebuilt casters with new components and wheels
- 1 year warranty (Domestic only)

Specifications

- "Brake not set" alarm
- Frame and mattress integration for superior therapy/safety performance
- FlexAfoot™ mechanism
- Enhanced turn assist (on A.I.R. surfaces only)
- Point of Care® controls - Large, easy-to-read design
- FlexAfoot™ mechanism
- Low position – surface edge to floor 18" with optional
- IntelliDrive® power transport 22"
- Maintains low height when placed in chair
- High position – surface center to floor 37"
- Headboard/push handles maintain same height regardless of patient surface height

Options

- A wide variety of room accessories
- IntelliDrive® power-assist transport system
- Digital head of bed angle indicator
- Advanced 3 mode bed exit
- Service Contracts for preventative maintenance
- 220-Volt transformers

