



# The DRE Integra AV-S Anesthesia System The Competition...



	IntegraAV-S(2 vap)	Fabius GS (2 vap)	Aestiva (3 vap)	Aespire (2 vap)
Dimensions" h x d x w	53"x 31" x 28"+abs	51"x 32" x 35"+abs	54"x 33" x 37"+abs	53"x 29" x 29"
Weight lb	220 including 2 vaps	260 including 2 vaps	387 including 3 vaps	268 including 2 vaps
Top shelf " w x d	28" x 19.5"	30" x 19"	26" x 16"	26" x 16"
Work surface " w x d	23" x 14"	27" x 13"	18" x 12"	23" x 14"
Writing tablet	Optional	yes pull out	Optional	Optional
Drawers" h x d x w	Up to 3 x (7.7" x 16" x 20")	3 x	1 x (4" x 10" x 15")	2 x (6.9" x 10" x 13")
Flow meter type	Dual flow tubes	LED & combined flow	Dual flow tubes	Optional Dual flow
Gases as standard	O2, N2O & Air	O2, N2O & Air	Air is optional	Air is optional
Cylinders as standard	1xO2,1xN2O	Optional	Optional	Optional
Electrical outlets x 4	Standard	Optional	Optional	Optional
O2 Therapy Flow meter	Optional	Optional	Optional	Optional
MRI compatible*	Standard*	No	No	No
Scavenging & Suction	Optional	Optional	Built in	Built in
Anti Hypoxic Device	Direct drive	S-ORC	Chain Link - 25	Chain Link - 25

Equipment for the way *you* operate

[www.dremed.com](http://www.dremed.com)