

Maquet 1130-1

Surgical Table

Equipment for the way *you* operate



Certified Refurbished

The C-Arm-compatible Maquet 1130-1 is ideal for radiography and intra-operative image intensification.

- Mobile operating table with a battery-powered, electrically-controlled hydraulic drive.
- Features an easy-to-use remote control and radio translucent table-top sections.
- Padded plates for upper, lower and seat sections are individually removable.
- Base and column cladding and the table frame are made of CrNi steel alloy that is resistant to disinfectants.
- Table frame has guide rails for X-ray cassettes.
- Rests on a T-shaped base with large, easy-rolling, ball-bearing castors for lateral and transverse motion. Castors are accessible for the side for disinfection and maintenance.



► Specifications

Dimensions: Length: 80";Width: 21"

Weight: 265 kg

Weight Limit: 500 lbs

Adjustment Ranges

- Trendelenburg: 30°
- Reverse Trendelenburg: 25°
- Lateral Tilt: 25° (right – left)
- Upper Back Section: 45° up, 65°down
- Lower Back: 65° up, 45°down
- Leg Plates: 26° up, 75°down

Dependable battery support: Maintenance-free batteries are integrated to the base of the Maquet 1130-1. The batteries can be used for two to three weeks of normal use (as many as 80 operations) without having to be charged.

Ref Code: 3011