



## PRODUCT FEATURES

- Twin pedal Trendelenburg
- Height adjustment clearly marked
- Pneumatic back adjustment
- Central locking casters
- Brake/steer pedals at each corner
- Fold down side rails
- Four infusion pole sockets
- Oxygen cylinder storage built-in
- 550 lb capacity
- 4" foam mattress
- Height adjustments from 23" to 26"
- Surface dimensions: 75" x 27" (l x w)
- Factory new with a 1-year warranty

## DRE Millennium 2 Hospital Stretcher

The DRE Millennium-2 Stretcher is a durable alternative for budget-conscious facilities.

Supporting up to 550-pounds on a four-inch foam mattress, the stretcher adjusts to Trendelenburg and Reverse-Trendelenburg positions. Infusion pole sockets at each corner eliminate the need for additional equipment in small spaces and ensure a secure transport. The back adjusts pneumatically for safety and ease of use.

Other standard features are central locking casters, brake/steer pedals at each corner, built-in oxygen cylinder storage, height adjustment (23 to 36 inches), and a one-year warranty.

