

Available from  **DRE**

Elegance · Economy · Ergonomics



ULTRASOUND IMAGING TABLES

- Bariatric strength
- Ergonomically advanced
- Customizable
- Built to last

Made in the USA 
with US & imported parts



MULTI-SPECIALTY

- **Ergonomic excellence** - all of the options in one table.
- **600 lb. Patient Weight Capacity** - dual lift column design ensures table stability and positioning flexibility



22" - 38" height range for convenient patient transfer and flexible working conditions



shown here with optional features

Powered fowler to easily lift your patient to a seated position or assist them in laying down.



PRODUCT SPECS

PATIENT WEIGHT CAPACITY	600 lbs. (272 kg.)
HEIGHT RANGE	22" - 38" (56-97 cm)
WIDTH	30" (76 cm)
LENGTH	74" (188 cm)
TOP	3 section; Fowler - 36" (91 cm), Seat - 24" (61 cm), Foot - 13" (33 cm)
TRENDELENBURG/REVERSE	±15°
PADDING	3" (8 cm) multi-layer
UPHOLSTERY	TerraTouch™ - over 20 colors
TABLE WEIGHT	175 lbs. (80kg.)
WARRANTY	3 year
CERTIFICATIONS	FDA listed, ETL approved, CE marked

PRODUCT FEATURES

X = INCLUDED O = OPTIONAL

CENTRAL BASE LOCKING SYSTEM	X
ELECTRIC FOWLER	X
ADJUSTABLE ACCESS/BACK SUPPORT (RIGHT)	X
ADJUSTABLE ACCESS/BACK SUPPORT (LEFT)	X
MANUAL ADJUSTABLE FOOT REST	X
HAND CONTROL (HEIGHT, TRENDELENBURG/REVERSE, FOWLER)	X
FOOT CONTROL (HEIGHT ONLY)	X
FOOT CONTROL (HEIGHT & FOWLER)	O
RETRACTABLE FOOT RESTS (STIRRUPS)	X
SIDE RAILS (2)	O
T-RAIL ADAPTER	O
DETACHABLE ARM BOARD	O
PAPER ROLL HOLDER	X
12" TABLE EXTENDER	O
QUICK LOCK HEAD REST & PAD	O
ACCESSORY MOUNTING PLATE (FOWLER END ONLY)	O



- **Elegance, Economy, Ergonomics**
a table built for the way you work.
- **Small base & table top design**
allows easy access to patient
- **Electric adjustments** allow you to bring your patient to you
- **Recessed side rails** for increased patient security, fold out of the way when not in use.



shown here with optional features



Adjustable back support with four positioning angles. Easily adjustable with one hand. Available on one or both sides.



Optional table extenders, such as the QuickLock™ Head Rest, allow accurate and comfortable cervical positioning.

PRODUCT FEATURES

X = INCLUDED O = OPTIONAL

CENTRAL BASE LOCKING SYSTEM	X
ELECTRIC FOWLER	X
ADJUSTABLE ACCESS/BACK SUPPORT (RIGHT)	X
ADJUSTABLE ACCESS/BACK SUPPORT (LEFT)	X
MANUAL ADJUSTABLE FOOT REST	X
HAND CONTROL (HEIGHT, TRENDELENBURG/REVERSE, FOWLER)	X
FOOT CONTROL (HEIGHT ONLY)	X
FOOT CONTROL (HEIGHT & FOWLER)	O
RETRACTABLE FOOT RESTS (STIRRUPS)	O
SIDE RAILS (2)	O
T-RAIL ADAPTER	O
PAPER ROLL HOLDER	X
12" TABLE EXTENDER	O
QUICK LOCK HEAD REST & PAD	O
ACCESSORY MOUNTING PLATE (FOWLER END ONLY)	O

PRODUCT SPECS

PATIENT WEIGHT CAPACITY	600 lbs. (272 kg.)
HEIGHT RANGE	22" - 38" (56-97 cm)
WIDTH	30" (76 cm)
LENGTH	74" (188 cm)
TOP	3 section; Fowler - 36" (91 cm) , Seat - 24" (61 cm), Foot - 13" (33 cm)
TRENDELENBURG/REVERSE	±15°
PADDING	3" (8 cm) multi-layer
UPHOLSTERY	TerraTouch™ - over 20 colors
TABLE WEIGHT	175 lbs. (80kg.)
WARRANTY	3 year
CERTIFICATIONS	FDA listed, ETL approved, CE marked





shown here with optional central base locking system



shown here with optional electric fowler



22" - 38" height range for convenient patient transfer and flexible working conditions

Optional T-Rail adapter can be attached at any of four locations on the table, allowing the use of any t-rail attachable accessory.



Optional table extender lengthens the table allowing for cervical access.



PRODUCT FEATURES

X = INCLUDED O = OPTIONAL

CENTRAL BASE LOCKING SYSTEM	O
ELECTRIC FOWLER	O
ADJUSTABLE ACCESS/BACK SUPPORT (RIGHT) AVAILABLE WITH FOWLER TOP ONLY.	O
ADJUSTABLE ACCESS/BACK SUPPORT (LEFT) AVAILABLE WITH FOWLER TOP ONLY.	O
HAND CONTROL (HEIGHT, TRENDLENBURG/REVERSE, FOWLER)	X
FOOT CONTROL (HEIGHT ONLY)	O
FOOT CONTROL (HEIGHT & FOWLER)	O
SIDE RAILS (2)	O
T-RAIL ADAPTER	O
PAPER ROLL HOLDER	X
12" TABLE EXTENDER	O
QUICK LOCK HEAD REST & PAD	O
ACCESSORY MOUNTING PLATE (FOWLER END ONLY)	O

PRODUCT SPECS

PATIENT WEIGHT CAPACITY	600 lbs. (272 kg.)
HEIGHT RANGE	22"- 38" (56-97 cm)
WIDTH	27" (69 cm) or 30" (76 cm)
LENGTH - FLAT TOP	Flat top - 72" (183 cm)
LENGTH-FOWLER TOP	Fowler- 34" (86 cm), Seat - 38" (97 cm)
TRENDLENBURG/REVERSE	± 15°
PADDING	3" (8 cm) multi-layer
UPHOLSTERY	TerraTouch™ - over 20 colors
TABLE WEIGHT	175 lbs. (80kg.)
WARRANTY	3 year
CERTIFICATIONS	FDA listed, ETL approved, CE marked

