

 Repotenciación Certificada

La solución neumática para cirugías en espacios pequeños

Las maquinas compactas Excel 110 y 210 vienen con O₂ y N₂O estándar. También cuenta con Sistema de Guarda Hipoxica y absorbedor GMS. La Excel 210 también incluye aire medicinal.



Foto:
Excel 210
Color y compartimientos pueden variar al de la foto.

Un vistazo rápido

Características

- » Capacidad para tres gases y dos vaporizadores
- » Ventilación de gran alcance y capacidad de monitoreo de lugares anestesiados secundarios

Opciones

- » Ventilador Electrónico 7800 con vaporizadores Tek 5
- » Ventilador 7000 con vaporizadores Tek 4
- » Vaporizador nuevo Sigma Delta

Specifications for the Excel 110

| | |
|-----------------------------|-------------------|
| Weight (approximate) | 270 lbs. (122 kg) |
| Options | |
| –Vaporizer | 17 lbs. (7.6 kg) |
| –Absorber | 21 lbs. (9.5 kg) |
| –Ventilator | 19 lbs. (8.5 kg) |

Dimensions

| | |
|--------|---|
| Height | 66" (168 cm) |
| | <i>from floor to upper most part of top shelf</i> |
| Depth | 30" (76 cm) |
| | <i>with casters and brake bar</i> |

Width

| | |
|----------------------------|-------------------|
| –Machine only | 20" (50.7 cm) |
| –Machine with monitor pod | 27" (68 cm) |
| –Absorber and Absorber Arm | Add 10" (25.5 cm) |

Shelving Systems

| | |
|---------------------|---|
| Middle Shelf | |
| –Maximum Load | 24 lbs. (11 kg) |
| –Usable height | 6" (15.2 cm) |
| –Shelf Size | 16.6" x 12.5" (42 cm x 31.7 cm) |
| | <i>varies slightly with shelf style</i> |

Top Shelf

| | |
|---------------|---|
| –Maximum Load | 48 lbs. (22 kg) |
| –Shelf Size | 21.5" x 14.2" (54.5 cm x 36 cm) |
| | <i>varies slightly with shelf style</i> |

Table-top

| | |
|-----------------------|---|
| –Height (above floor) | 33.5" (84 cm) |
| –Size | 20" x 11.7" (50.8 cm x 29 cm) |
| | <i>varies slightly with shelf style</i> |

Drawer Cabinet

- Contains two 15.5" x 10.5" x 4" (39.4 cm x 26.7 cm x 10.2 cm) drawers

Absorber Post Mounting Assembly

| | |
|-----------------------------|---------------------------------------|
| –Absorber Arm Length | 14" (35.5 cm) |
| | <i>varies slightly with arm style</i> |
| –Vertical Adj (above floor) | 10" to 26" (25.4 cm to 66 cm) |
| –Horizontal Adj | 6.5" to 11" (16.5 cm to 27.9 cm) |

Casters

| | |
|----------|--|
| Diameter | 5" (12.7 cm) |
| | <i>front casters have a foot-operated brake-bar lock</i> |

Common Gas Outlet

- Equipped with a latching, positive engagement, bayonet type connector. The common gas outlet connector will also accept standard 22 mm OD or 15 mm ID conical friction fit connectors.

Specifications for the Excel 210

| | |
|-----------------------------|-------------------|
| Weight (approximate) | 300 lbs. (136 kg) |
| Options | |
| –Vaporizer | 17 lbs. (7.6 kg) |
| –Absorber | 21 lbs. (9.5 kg) |
| –Ventilator | 19 lbs. (8.5 kg) |

Dimensions

| | |
|--------|---|
| Height | 66" (168 cm) |
| | <i>from floor to upper most part of top shelf</i> |
| Depth | 30" (76 cm) |
| | <i>with casters and brake bar</i> |

Width

| | |
|----------------------------|-------------------|
| –Machine only | 22" (55.8 cm) |
| –Machine with monitor pod | 29" (73.6 cm) |
| –Absorber and Absorber Arm | Add 10" (25.5 cm) |

Shelving Systems

| | |
|---------------------|---|
| Middle Shelf | |
| –Maximum Load | 24 lbs. (11 kg) |
| –Usable height | 6" (15.2 cm) |
| –Shelf Size | 18.6" x 12.5" (47.2 cm x 31.7 cm) |
| | <i>varies slightly with shelf style</i> |

Top Shelf

| | |
|---------------|---|
| –Maximum Load | 48 lbs. (22 kg) |
| –Shelf Size | 23.5" x 14.2" (59.6 cm x 36 cm) |
| | <i>varies slightly with shelf style</i> |

Table-top

| | |
|-----------------------|---|
| –Height (above floor) | 33.5" (84 cm) |
| –Size | 21.7" x 11.7" (55 cm x 29 cm) |
| | <i>varies slightly with shelf style</i> |

Drawer Cabinet

- Contains two 15.5" x 10.5" x 4" (39.4 cm x 26.7 cm x 10.2 cm) drawers

Absorber Post Mounting Assembly

| | |
|-----------------------------|---------------------------------------|
| –Absorber Arm Length | 14" (35.5 cm) |
| | <i>varies slightly with arm style</i> |
| –Vertical Adj (above floor) | 10" to 26" (25.4 cm to 66 cm) |
| –Horizontal Adj | 6.5" to 11" (16.5 cm to 27.9 cm) |

Casters

| | |
|----------|--|
| Diameter | 5" (12.7 cm) |
| | <i>front casters have a foot-operated brake-bar lock</i> |

Common Gas Outlet

- Equipped with a latching, positive engagement, bayonet type connector. The common gas outlet connector will also accept standard 22 mm OD or 15 mm ID conical friction fit connectors.