



Certified Refurbished



Your pneumatic solution for surgeries in **close quarters**

The compact Excel 110 and 210 come standard with O₂ and N₂O, as well as a Hypoxic Guard and a GMS absorber. The Excel 210 also includes medical air.

Pictured: Excel 210
Color and drawer pack configurations may vary from picture.

► A Quick Look

Features

- Accommodates three gases and two vaporizers
- Powerful ventilation and monitoring capabilities for secondary anesthetizing locations.

Options

- 7800 Electronic Ventilator with Tek 5 Vaporizers
- 7000 Ventilator with Tek 4 Vaporizers
- New Sigma Delta Vaporizer

Ohmeda Excel 110/210

Anesthesia Machine

Equipment for the way *you* operate

Specifications for the Excel 110

Weight (approximate)	270 lbs. (122 kg)
Options	
–Vaporizer	17 lbs. (7.6 kg)
–Absorber	21 lbs. (9.5 kg)
–Ventilator	19 lbs. (8.5 kg)

Dimensions

Height	66" (168 cm)
	<i>from floor to upper most part of top shelf</i>
Depth	30" (76 cm)
	<i>with casters and brake bar</i>

Width

–Machine only	20" (50.7 cm)
–Machine with monitor pod	27" (68 cm)
–Absorber and Absorber Arm	Add 10" (25.5 cm)

Shelving Systems

Middle Shelf	
–Maximum Load	24 lbs. (11 kg)
–Usable height	6" (15.2 cm)
–Shelf Size	16.6" x 12.5" (42 cm x 31.7 cm)
	<i>varies slightly with shelf style</i>

Top Shelf

–Maximum Load	48 lbs. (22 kg)
–Shelf Size	21.5" x 14.2" (54.5 cm x 36 cm)
	<i>varies slightly with shelf style</i>

Table-top

–Height (above floor)	33.5" (84 cm)
–Size	20" x 11.7" (50.8 cm x 29 cm)
	<i>varies slightly with shelf style</i>

Drawer Cabinet

- Contains two 15.5" x 10.5" x 4" (39.4 cm x 26.7 cm x 10.2 cm) drawers

Absorber Post Mounting Assembly

–Absorber Arm Length	14" (35.5 cm)
	<i>varies slightly with arm style</i>
–Vertical Adj (above floor)	10" to 26" (25.4 cm to 66 cm)
–Horizontal Adj	6.5" to 11" (16.5 cm to 27.9 cm)

Casters

Diameter	5" (12.7 cm)
	<i>front casters have a foot-operated brake-bar lock</i>

Common Gas Outlet

- Equipped with a latching, positive engagement, bayonet type connector. The common gas outlet connector will also accept standard 22 mm OD or 15 mm ID conical friction fit connectors.

Specifications for the Excel 210

Weight (approximate)	300 lbs. (136 kg)
Options	
–Vaporizer	17 lbs. (7.6 kg)
–Absorber	21 lbs. (9.5 kg)
–Ventilator	19 lbs. (8.5 kg)

Dimensions

Height	66" (168 cm)
	<i>from floor to upper most part of top shelf</i>
Depth	30" (76 cm)
	<i>with casters and brake bar</i>

Width

–Machine only	22" (55.8 cm)
–Machine with monitor pod	29" (73.6 cm)
–Absorber and Absorber Arm	Add 10" (25.5 cm)

Shelving Systems

Middle Shelf	
–Maximum Load	24 lbs. (11 kg)
–Usable height	6" (15.2 cm)
–Shelf Size	18.6" x 12.5" (47.2 cm x 31.7 cm)
	<i>varies slightly with shelf style</i>

Top Shelf

–Maximum Load	48 lbs. (22 kg)
–Shelf Size	23.5" x 14.2" (59.6 cm x 36 cm)
	<i>varies slightly with shelf style</i>

Table-top

–Height (above floor)	33.5" (84 cm)
–Size	21.7" x 11.7" (55 cm x 29 cm)
	<i>varies slightly with shelf style</i>

Drawer Cabinet

- Contains two 15.5" x 10.5" x 4" (39.4 cm x 26.7 cm x 10.2 cm) drawers

Absorber Post Mounting Assembly

–Absorber Arm Length	14" (35.5 cm)
	<i>varies slightly with arm style</i>
–Vertical Adj (above floor)	10" to 26" (25.4 cm to 66 cm)
–Horizontal Adj	6.5" to 11" (16.5 cm to 27.9 cm)

Casters

Diameter	5" (12.7 cm)
	<i>front casters have a foot-operated brake-bar lock</i>

Common Gas Outlet

- Equipped with a latching, positive engagement, bayonet type connector. The common gas outlet connector will also accept standard 22 mm OD or 15 mm ID conical friction fit connectors.