



We Have Other Models!

Pedigo 2030

BLANKET & FLUID WARMER

PRODUCT FEATURES

- New and Fully insulated
- Compartment size is: 7.7 cubic feet -or- 18-7/8" w x 26-1/2" d x 26-3/4" h
- Overall size is: 24-1/4" w x 33-3/4" d x 41-3/4" h
- Easy inventory visualization through optional glass doors
- Electrothermal radiant heat
- Safety coated interior
- Individual controls adjust to 200°F
- 125 VAC or 230 VAC
- UL & CSA approved
- Casters or fixed legs
- Can be ordered in a stackable design with another P-2020 or with a P-2120 Fluid Warming Cabinet

AN UPGRADE TO THE DELUXE MODEL OFFERS:

4-Digit LED Display, (Additional digit allows for 24-hour time programming and Fahrenheit/Celsius display), Control Lock Programming, (Ability to lock temperature into unit. Auto Time Programming (Allows the user to program the controls to automatically turn on and off once during a 24-hour period at selected times.) Prevents inadvertent changing of the temperature setting.), One-Touch Temperature Increase/Decrease buttons, Rounded Front Face Panel, and Offset Door Handle.

OTHER MODELS & SIZES

- P-2010 16-3/16" w x 25" d x 10-5/6" h
- P-2040 24-1/4" w x 33-3/4" d x 74 " h
- P-2050 39" w x 30" d x 67-1/4" h

