



REFURBISHED
SKYTRON 6001
 Surgical Table

- General purpose surgical table with remote control and radio tops
- C-arm interface capabilities meet all imaging requirements
- The stability of the table and rigidity of the top provide assurance that the patient position will be maintained as desired
- Optional foot control

SPECIFICATIONS:

- | | |
|-------------------|---------|
| • Weight Capacity | 500 lbs |
| • Min. Height | 28" |
| • Max. Height | 45" |
| • Table Width | 20" |
| • Table Length | 75" |
| • Lateral Tilt | 30° |
| • Trendelenburg | 30° |
| • Top Rotation | 180° |