



Certified Refurbished



The Skytron 6500 is C-Arm compatible and has a 500-lb. patient capacity

- Manual, electric, and remote control models available
- C-Arm compatible and X-ray top set
- 180° table top rotation
- Removable powered leg section
- Positions include Trendelenburg, back (raise or lower), side (lateral tilt), foot (raise or lower), and flex ("v" or inverted "v")
- Optional accessories include armboards, stirrups, anesthesia screens, and x-ray top
- Accommodates taller patients with the adjustable head section
- Use the reversible head section as a foot extension while transferring patients
- Operates on line or battery power
- Sturdy, easy-to-clean stainless steel column and base
- Brakes are completely mechanical and provide dependable performance, even when there is no power to the machine
- Helps reduce surgeon back strain via a narrow base design that has ample toe space clearance
- Push button break system
- Large perineal cut-out
- Kidney bridge with 5" of additional lift

► Specifications

- Weight Lifting Capacity: 500 lbs.
- Maximum Height: 47"
- Minimum Height: 26"
- Reflex: 90°
- Flex: 70°
- Foot / Leg Down: 105°
- Foot / Leg Up: 20°
- Back Down: 40°
- Back Up: 90°
- Trendelenburg: 30°
- Reverse Trendelenburg: 30°
- Tilt Right: 30°
- Tilt Left: 30°
- Power Requirements: 120 VAC, 60 HZ 300 Watts
- Current Leakage: Less than 100 Micro Amps
- Power Cord: 15 Feet with Hospital Grade Connector