

Skytron 6500-HD

Surgical Table

Equipment for the way *you* operate



Certified Refurbished



Capable of lifting and articulating a full 850-lb. patient load

► Specifications

Physical Specifications

- Weight Capacity: 850 lbs.
- Height Adjustment : 27" to 44"
- Top Rotation: 180°
- Foot / Leg Down: 105° (Removable)
- Foot / Leg Up: 20°
- Back Down: 40°
- Back Up: 60°
- Trendelenburg and Reverse Trendelenburg: 20°
- Lateral Tilt: 20°

Electrical Specifications

- Power Requirements: 120 VAC, 60 HZ 300 Watts
- Current Leakage: Less than 100 Microamps
- Power Cord: 15 Feet with Standard Duplex Connector

- General use electro-hydraulic surgical table
- Can lift and fully articulate patients weighing up to 850 lbs.
- C-Arm compatible and X-ray top set
- 180° table top rotation
- Buttons are clearly labeled, making the table easy to use and enhancing patient safety
- Fail-safe feature makes sure brakes are locked prior to table use
- Stainless steel cover is attractive and easy to clean
- Self-leveling, four-point brake system
- Built-in kidney bridge
- Manually-operated head section features a full width quick-release bar
- Narrow base