

New, out of the box Premier 4500 stretchers

**Order now for special discount pricing**

(limited quantities available)



**Includes a Four-Year Manufacturer's Warranty**

### Specifications

- Height Adjustment: 22" to 35.5" (56 cm to 90 cm)
- Integrated Fold Down/Collapsible Side Rails: 60.6" x 12.5" (154 cm x 32 cm)
- Raise/Lower Pedals: Side Access
- Hands-Free Trendelenberg/Reverse Trendelenberg: +18° to -18°
- 8" (20 cm) All Surface, Four-Wheel Locking Castors
- Accu-Steer Fifth Wheel Steering
- Techlem Posi-Lock Braking System
- Pneumatic Assist Head Section: 0° to 90°
- Wraparound Bumper
- Weight Capacity: 500 lbs.

### Dimensions

- Total Dimensions: 81" x 32.5" (206 cm x 83 cm)
- Litter Dimensions: 81" x 28.5" (206 cm x 73 cm)
- Mattress Dimensions: 75" x 28.5" x 3" (191 cm x 73 cm x 7.5 cm)

